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CANADA Mölnlycke Health Care Inc.,
2010 Winston Park Drive, Suite 100,
Oakville, Ontario, L6H 5R7
Tel: +1 905 829 1502

Issued 2019-04

Master 499017 Rev. 08

327322-01

PD-556430 Rev. 00



Made in UK

Mölnlycke®
Exufiber® Ag+

Gelling Fibre Dressing with Silver

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www.molnlycke.com

Exufiber® Ag+

Gelling Fibre Dressing with Silver

Rx only Federal [US] Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

Product description

Exufiber® Ag+ wound dressing is a sterile, soft, nonwoven pad or ribbon. Exufiber Ag+ consists of highly absorbent and gel forming polyvinyl alcohol (PVA) nonwoven fibres. The nonwoven pad or ribbon is coated with silver sulphate on both sides. The silver is released within the dressing in contact with fluid, and acts as a preservative in the dressing to inhibit or reduce microbial growth in the dressing.

In vitro testing demonstrates that Exufiber Ag+ kills a wide range of known wound pathogens and creates an antimicrobial environment within the dressing. Exufiber Ag+ gives a x4 log reduction, in the dressing, of the following Gram positive bacteria, Gram negative bacteria and fungi through seven (7) days in vitro studies. By reducing the number of micro-organisms, Exufiber Ag+ may also reduce dressing odour.

	7 days sustained activity test (simulated wound model with serum)
Group	Strain
Gram positive species	<i>Enterococcus faecalis</i> [VRE]
	<i>Staphylococcus aureus</i> [MRSA]
	<i>Staphylococcus epidermidis</i>
	<i>Streptococcus pyogenes</i>
Gram negative species	<i>Acinetobacter baumannii</i>
	<i>Klebsiella pneumoniae</i>
	<i>Enterobacter cloacae</i>
	<i>Pseudomonas aeruginosa</i>
Fungal species	<i>Candida guilliermondii</i>
	<i>Candida tropicalis</i>
	<i>Candida albicans</i>
	<i>Aspergillus fumigatus</i>

Exufiber Ag+ has shown to have an unchanged microbial reduction through seven (7) days when tested in vitro with a secondary dressing.

Exufiber Ag+ maintains a moist wound environment in combination with a secondary dressing.

Exufiber Ag+ can be used under compression bandaging.

Indications

Exufiber Ag+ is intended to be used on the following medium to high exuding wounds:

- Leg ulcers (venous stasis ulcers, arterial ulcers and ulcers of mixed etiology) and diabetic foot ulcers
- Pressure ulcers (partial and full thickness)
- Partial thickness burns (second degree)
- Donor sites and other wounds that are prone to bleeding, such as debrided wounds (mechanically or surgically)
- Traumatic wounds
- Surgical wounds that heal by primary intention, such as dermatological and surgical incisions (e.g. orthopaedic and vascular), and surgical wounds left to heal by secondary intention, such as dehisced surgical incisions
- Oncology wounds with exudate, such as fungoides-cutaneous tumours, fungating carcinoma, cutaneous metastasis, Kaposi’s sarcoma and angiosarcoma

Exufiber Ag+ may be used for management of wounds as an effective barrier to bacterial penetration of the dressing, as this may help to reduce the risk of infection.

Indicated wear time: up to seven (7) days.

Instructions for use

Shallow wounds

Application:

Cleanse the wound with saline solution or water according to clinical practice.

- Dry the surrounding skin thoroughly.
- Choose the correct size of dressing to be able to cover the entire wound. For best result, Exufiber Ag+ should overlap the dry surrounding skin by at least 1-2cm for the smaller sizes (up to 10x10cm) and 5cm for larger sizes. The dressing will swell, causing contraction of the dressing edges, as it absorbs wound fluid and starts gelling.
- Apply a dry Exufiber Ag+ dressing to the wound.
- Fixate with an appropriate secondary dressing. The choice depends on the exudate level.
- Compression therapy may be used in conjunction with Exufiber Ag.

Deep wounds and cavities

Application:

Cleanse the wound with saline solution or water according to clinical practice.

- Dry the surrounding skin thoroughly.
- Loosely pack ribbon or pad into the wound to allow room for swelling of the dressing. When using the ribbon dressing, cut to appropriate length leaving a small overhang of 2-3cm outside the wound for easy retrieval.
- Cover with an appropriate secondary dressing. The choice depends on the exudate level.
- Fixate (or cover) with an adhesive dressing or a bandage as appropriate.

Donor sites

Application:

Harvest tissue and achieve hemostasis (according to normal routine).

- Choose the correct size of sheet dressing to be able to cover the entire wound. The Exufiber Ag+ dressing should overlap the dry surrounding skin by at least 5cm, the dressing will swell, causing contraction of the dressing edges, as it absorbs wound fluid and starts gelling.
- Apply a dry Exufiber Ag+ dressing to the wound.
- Cover with an appropriate secondary dressing. The choice depends on the exudate level.
- Fixate (or cover) with an adhesive dressing or a bandage as appropriate.

Dressing change and removal

Change the Exufiber Ag+ dressing when saturation is reached. Exufiber Ag+ can be left in place for up to 7 days, depending on wound condition or as indicated by clinical practice.

- Remove the secondary dressing and discard in appropriate way.
- Remove the Exufiber Ag+ dressing by gently cleansing/flushing with saline solution or water according to clinical practice. Any non-gelled material will moisten in contact with the saline.
- If the dressing dries out and is difficult to remove, it should be moistened with saline solution or water according to local policies and allowed to soak until it lifts easily. It may take several minutes for Exufiber Ag+ to transform into a gel. Remove the dressing by gently cleansing/flushing.

Contraindications

- Exufiber Ag+ is not intended for dry wounds or full thickness burns.
- For external use only.

Precautions

- All wounds should be inspected frequently. In case of signs of clinical infection, consult a health care professional. Exufiber Ag+ does not replace the need for systemic therapy or other adequate infection treatment.
- Exufiber Ag+ should be used under the supervision of a qualified health care professional.
- Do not use on patients with a known sensitivity to silver or any other contents of the dressing.
- Clinicians/ Healthcare Professionals should be aware that there are very limited data on prolonged and repeated use of silver containing dressings, particularly in children and neonates.
- Frequent or prolonged use of this product may result in permanent discoloration of skin.
- Prior to commencing radiation therapy remove Exufiber Ag+ if product is present in the treatment area. A new dressing can be applied following treatment.
- Avoid contact with electrodes or conductive gels during electronic measurements, e.g. electrocardiograms (ECG) and electroencephalograms (EEG).
- Do not use Exufiber Ag+ together with oxidising agents such as hypochlorite solutions or hydrogen peroxide.
- Other than saline solution or water, the interaction of cleansing agents in combination with Exufiber Ag+ has not been demonstrated.
- The interaction of Exufiber Ag+ with topical treatments has not been demonstrated.

Storage and disposal

Exufiber Ag+ should be stored in dry conditions below 25°C/77°F and protected from direct sunlight. Note that any colour variances in Exufiber Ag+ do not affect product performance or safety.

Disposal should be handled according to local environmental procedures.

Other information

Exufiber Ag+ contains 0.2 mg/cm² silver, this translates to the below amounts in the different dressing sizes:

Dressing sizes		Total silver content
cm	inches	mg
5x5	2x2	5
10x10	4x4	20
10x12	4x4.8	24
15x15	6x6	45
20x30	8x12	120
2x45 (ribbon)	0.8x17.7 (ribbon)	18

The PVA fibres used in the product may change colour when it is exposed to light, air and/ or heat. The colour change has no influence on product properties when used before expiry date.

Do not reuse. If reused performance of the product may deteriorate, cross contamination may occur.

Sterile. Do not use if inner package is damaged or opened prior to use. Do not re-sterilise.

Exufiber is a registered trademark of Mölnlycke Health Care AB.

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